

NOVAPED PROTECT ORTHOPAEDIC FOOT SUPPORTS FOR WORK SAFETY SHOES



ELTEN

atlas
SAFETY SHOES

STEITZ SECURA
SICHERHEITSSCHUHE

HAIX

Holík

schein orthopädie
service
— Germany since 1879 —

Wehelp

INCREASING QUALITY OF LIFE
AND WELLBEING OF PEOPLE!

NOVAPED PROTECT ORTHOPAEDIC FOOT SUPPORTS FOR WORK SAFETY SHOES

According to DGUV regulations 112-191, only orthopaedic insoles with a valid type approval may be put into safety shoes for them to continue to comply with the EN 20345 standard.

NovaPED protect orthopaedic foot supports have been specially developed for use in safety shoes from the specified safety shoe manufacturers. Thanks to the use of certified cover materials and compliance with the corresponding manufacturing instructions, the NovaPED protect orthopaedic foot supports are compliant with the type approval certifications of the safety shoe manufacturers. A list of certified safety shoes can be requested or downloaded from the safety shoe manufacturers.

Safety and quality are combined in the NovaPED protect insoles. High-quality materials that meet the occupational safety requirements and the requirements of the MDR and are processed to high standards are used. The NovaPED protect orthopaedic foot supports can be further processed quickly and easily, in strict compliance with the manufacturing instructions in order to maintain conformity with the type approval certificates.



Scan the QR code to conveniently download the manufacturing instructions or go straight to WWW.SCHEIN-PROTECT.DE


NOVAPED PROTECT S90

Supporting insoles produced using the film technique



NOVAPED PROTECT SOFT

Bedding orthopaedic foot supports produced using soft cushioning technique



NOVAPED PROTECT S90 **ELTEN**

Supporting insoles produced
using the film technique



Certified cover material
ELTEN ESD PRO

Art. no.	Description	Size
989610000	Novaped Protect S90 ELTEN	35-48
989610800	Novaped Protect S90 ELTEN, loose cover	35-48

NOVAPED PROTECT SOFT **ELTEN**

Bedding orthopaedic foot supports
produced using soft cushioning technique



Certified cover material
ELTEN ESD PRO

Art. no.	Description	Size
979610000	Novaped Protect soft ELTEN	35-50
979610011	Novaped Protect soft ELTEN, partially adhered	35-50

COVER MATERIAL **ELTEN**

Art. no.	Description	Size
011961	ELTEN ESD PRO certified cover material, 4.3 mm	36-38 39-40 41-42 43-44 45-46 47-52



NOVAPED PROTECT S90 *atlas*

Supporting insoles produced using the film technique



Certified cover material
ATLAS® Ergo-Med

Art. no.	Description	Size
959614000	Novaped Protect S90 ATLAS®	35–48
959614550	Novaped Protect S90 ATLAS®, loose cover	35–48

NOVAPED PROTECT SOFT *atlas*

Bedding orthopaedic foot supports produced using soft cushioning technique



Certified cover material
ATLAS® Ergo-Med

Art. no.	Description	Size
979614000	Novaped Protect soft ATLAS®	35–50
979614011	Novaped Protect soft ATLAS®, partially adhered	35–50

COVER MATERIAL *atlas*

Art. no.	Description	Size
01.1964000	ATLAS® Ergo-Med certified cover material	335 x 240 x 2,5 mm




NOVAPED PROTECT SOFT STEITZ SECURA

Bedding orthopaedic foot supports produced using soft cushioning technique



Certified cover material
STEITZ SECURA ERGO SOFT ESD with forefoot perforation

Art. no.	Description	Size
979612000	Novaped Protect soft STEITZ SECURA	35-50
979612011	Novaped Protect soft STEITZ SECURA, partially adhered	35-50
011962000	STEITZ SECURA ERGO SOFT ESD certified cover material	390 x 290 x 2 mm



NOVAPED PROTECT SOFT HAIX

Bedding orthopaedic foot supports produced using soft cushioning technique



Certified cover material
HAIX CO-System

Art. no.	Description	Size
979613000	Novaped Protect soft HAIX	35-50
979613011	Novaped Protect soft HAIX, partially adhered	35-50
011963000	HAIX CO-System certified cover material	360 x 245 x 3,7 mm



NovaPED PROTECT

by **schein**

The other NovaPED protect orthopaedic foot supports, which are type-tested in conjunction with various safety shoes by an independent testing institute in accordance with EN 20345, also ensure increased made-to-measure safety.

They meet all the required test criteria and their contours are adapted to the safety shoes.

NOVAPED PROTECT SOFT Holik

Bedding orthopaedic foot supports produced using soft cushioning technique



Orthopaedic foot support, certified for Holik firefighting boots

Art. no.	Description	Size
979611000	Novaped Protect soft Holik	35-50
979611011	Novaped Protect soft Holik, partially adhered	35-50

NovaPED work

by **schein**

The NovaPED work insole range offers further options in terms of orthopaedic foot supports for additional safety shoes.



Schein Orthopädie Service KG
 Hildegardstr. 5
 42897 Remscheid
 Germany
 Tel. +49 2191 910-200
 Fax +49 2191 910-100
 service@schein.de
 www.schein.de



099993216

