


SCHEINWORKS morph






SCHEINWORKS morph

Your perfect craftsmanship – digitally prepared by us

Scheinworks morph newly developed piece of software that allows you to digitally design your individual orthopaedic foot supports in the fastest way possible. What do you need to do? Upload 2D or 3D scan data of your patient's feet, enter the patient's complaints and choose the appropriate product.

Our automatic morphing process combines all the information and suggests an individual orthopaedic foot support to you on that basis. This suggestion is already adapted to the shape of the foot, the indication and the desired correction strength. The result is a perfect and 100% individual insole surface!


Of course, you are then free to adjust and change the insole according to your wishes and preferences. Finally, you can choose whether the orthopaedic foot support should be printed or milled – and that's it! We receive the order and produce it immediately




Natural shapes – perfectly supported by technology

The benefits

- Reduced processing time due to individual and problem-specific insole suggestion from the software, based on the indication
- Individual post-processing of the digital insole suggestion
- Option to use print images, 2D scans or 3D scans
- Intuitive to use
- View of the morphed insole in 3D
- Particularly customisable design
- Order history and customer database (overview/status, etc.)
- No warehouse storage



From the foot straight
to the individually
morphed insole



04.0949.03.06.22



099993409



Wehelp

WHATSAPP:

+49 2191 910-200



Schein Orthopädie Service KG
Hildegardstr. 5
42897 Remscheid
Germany
Tel. +49 2191 910-200
Fax +49 2191 910-100
service@schein.de
www.schein.de